



VARIOUS AIR
OUTLET MODULES
MAX
-20°C

MACURAY

MACURAY
2 in 1 Air Cooling & Evacuator

MACURAY

2 in 1 System, Cryo air cooler & Evacuator system



Cooling Air

Smoke Evacuator

MACURAY


Two in One MACURAY
2 functions can be used simultaneously

Macuray maintains the low temperature with cooling air to relieve pain. In addition, Macuray has a smoke evacuator function which designed to effectively absorb and filter surgical smoke, odors, particulates, and other potentially dangerous vapors during face and body treatment.

2 in 1 System Compatible With Various Skin Treatments


Cooling Air

Various Air Outlet Modules
MAX -20 °C (±5 °C)



Smoke Evacuator

Easy On/Off Control by Foot Switch
Adopt a High performance Smoke Filter



Specifications

Power	220~230V AC 50/60Hz	Range of Temperature	Cooling Air : Max -20 °C (±5 °C) (Based on 23°C Surrounding Temp & Humidity is 40%)
Current Consumption	Cooling Device : 450W Smoke Evacuator : 1,100W	Dimension	340(W) x 400(D) x 910(H)mm, 1310(H)mm with Pole (Available separation Pole)
Cooling Air Delivery	45 CFM	Weight	35 kg (Net.)
Vacuum Power	2,400 mmH2O		