

[Q&A]

Q Why is Xyblo Special?

- A**
- Xyblo is Korea's First, Esthetic First in incorporating Quad-Polar technology.
 - It uses 0.5 MHz, 1 MHz cross-modulated high frequency.
 - It features a cross output of high frequency and HVMC (High Voltage Multi Current).
 - It is possible to care for the epidermis, dermis, subcutaneous fat, and fascia.

Q Is the pain severe?

- A** Xyblo is the most painless procedure among lifting procedures.

Q What is the treatment cycle?

- A** X-mode is recommended 10X once or twice a week
DPL is recommended 10X once a week.

Q Is there anything special about the procedure?

- A** Xyblo is complete with exosomes. When combined with exosomes, you can see more excellent effects.



Xyblo

💡 Recommend it to these people!

- People who want management of their wrinkles near the eye and mouth
- People who want management of their skin texture and skin tone
- People who are worrying over their sagging, non-elastic skin
- People who are worrying because of their double chin
- People who want management of their face contour



Xyblo

IN and OUT at once,
the secret is X-mode!
3D CONTOURING

Xyblo

Wrinkle Management / Skin Texture, Skin Tone /
Elasticity Management / Double Chin




The Principles of Xyblo


What is Xyblo?

With celebrity lifting, it simultaneously implements painless reduced lifting and full lifting!

- Esthetic's First Implementation of Quad-Polar Technology
- 0.5MHz, 1MHz Cross-Modulated High Frequency HVMC(High Voltage Multi Current) Electrical Stimulation



- 1 HIFU : SMAS(Fascia Layer) Stimulation
- 2 0.5MHz RF : 2~3.5mm
- 3 1.0MHz RF : 0.2~2.5mm
- 4 HVMC : 0.2~5mm
- 5 Pain-free multi-skin regeneration
- 6 Epidermis, dermis, subcutaneous fat, fascia regeneration effect
- 7 Collagen regeneration, improvement of wrinkles and sagging
- 8 Minimize side effects such as bruising and erythema during the procedure
- 9 Down Time minimization



Xyblo's Characteristics

Safe Procedure

Esthetic's First Implementation of Quad-Polar Technology

Multi Skin Rejuvenation

Skin Texture, Elasticity, Sagging Skin, Double Chin Improvement

Sophisticated Operation

- 1 Double Chin and Contour with HIFU
- 2 Face Elasticity and Collagen Production with HVMC
- 3 Detailed treatment is possible with HIFU handpiece and Quad-Polar handpiece

Elasticity Rising From the Skin

By delivering thermal RF energy to the dermis layer, it stimulates the dermal layer tissue to activate collagen and elastin regeneration, making the skin firm and transparent.

- 1 Cross-modulation System of RF and HVMC
- 2 0.5MHz and 1MHz Cross-Modulation System
- 3 Cross-Modulation System of Exo-Boost and ESW (Electric Shock Wave)
- 4 Separate mode available for Exo-Boost
- 5 All functions available separately
- 6 Equipped with Auto Function that is easy for the operator to use

Xyblo



INDICATION

Wrinkles Near Eye and Mouth

Skin Texture, Skin Tone

Skin Elasticity

Double Chin

Face Contour

Sagging Skin

Body Line

IN and OUT at once, the secret is X-mode!
3D CONTOURING
It's Possible Because It's Xyblo!