


1Q CATHETER

INTRODUCTION (EN)

1Q CATHETER


1Q SAFETY GUIDED CATHETER



Filler Cannula&Puncture Needle

Advantages

- It is easy to find the treatment area.
- The filler cannula can easily be inserted into the punching needle, which minimizes skin damage.
- There is almost no difference between the outer diameter of the filler cannula and the inner diameter of the punching needle, which increases the precision of the procedure.
- It reduces treatment time and treatment fatigue.



Procedure



1. Use the enclosed punching needle to perforate the treatment area.



2. Insert the filler cannula into the punching needle.



3. Insert the filler cannula into the punching needle, and inject the filler as it is.

Specification

Packing unit: 50pcs/ BOX



Filler Cannula		Puncture Needle	
Gauge	Length	Gauge	Length
22G	50mm	21G	5mm
	70mm		
23G	50mm	22G	5mm
	70mm		
25G	40mm	23G	5mm
	50mm		
	70mm		
27G	40mm	25G	5mm
	50mm		
	70mm		
30G	25mm	25G	5mm





Tel : +971 6 55 79 818
info@maria-trading.com
www.maria-trading.com
Sharjah - UAE

