

Dr. Cannula

- **Dr. Cannula**



- It is safer than the regular needle because there are holes on the side of the cannula needle then you can do more uniform injection.


- Got FDA Thailand

- Dr. Cannula



- **Marked scales** on the cannula needle. So The physician performs the procedure at the correct region of treatment.
- There are **holes on the side of the Cannula** needle can inject filler wide and uniformly across the side.
- It can be injected into a wide area with only one insertion because cannula **needle is flexible.**

- **Dr. Cannula**



DR.Cannula Normal needle

<Dr. Cannula>

It is inserted smoothly to minimize nerve damage and swelling because the needle tip is round (Blunt).
In addition, it can be performed without causing any damage to the blood vessels and nerves.

- **Dr. Cannula**



DR.Cannula Normal needle

<Dr. Cannula>

Screw type hub :Accurately and tightly connecting with syringe.

<Third-party products>

Normal type hub :The syringe and needle are easily separated in case of injecting the high pressure or the high density solution injection so that the contents are leaking.

- **Dr. Cannula**

Needle Gauge	Needle Length		Outer Diameter	Inner Diameter
21G	50mm		0.80mm	0.58mm
23G	50mm		0.60mm	0.42mm
25G	40mm	50mm	0.50mm	0.32mm
27G	40mm	50mm	0.40mm	0.27mm

-Specification and Packing unit : 40pcs / 1BOX, Individually sterilized packaging

Thank you for reading



Tel : +971 6 55 79 818
info@maria-trading.com
www.maria-trading.com
Sharjah - UAE

